


**VII Curso de Actualización en Oncopatología y III Curso-Taller de Citometría de Flujo**  
**Facultad de Medicina de la Universidad de Concepción**  
**02 y 03 de mayo del 2019, Concepción, Chile**


# **Enfermedad mínima residual en Mieloma Múltiple**

**Evangelina Agriello**  
**Bahía Blanca, Argentina**

# Plasma cell dyscrasias (PCD)




# Plasma cell dyscrasias (PCD)



# Laboratory SOPs

**Validated by an independent consortium of clinical flow cytometry experts – EuroFlow™**

- Sample collection and transport
- Instrument set up, calibration and compensation
- Sample processing and staining







## PRE-ANALYTICAL PHASE

- SPECIMEN HANDLING
- ANTICOAGULANT
- COLLECTION TIMING
- PRESERVATION  
CONDITIONS

## ANALYTICAL PHASE

- FLUOROCHROME SELECTION
- ANTIBODY COMBINATIONS
- SAMPLE PREPARATION
- INSTRUMENT SET-UP
- RAW DATA GENERATION

## POST-ANALYTICAL PHASE

- DATA INTEGRATION
- REPORTING

- Panel design
  - Fluorochromes
  - Antibodies
- Sample preparation and staining
  - Staining of surface markers only
  - Staining of surface markers including Igs
  - Combined staining of Intracellular and Surface membrane markers
  - Low cellularity samples/rare populations – The Bulk lyse protocol
- SOPs for instrument settings
  - FSC, SSC
  - Fluorescence channels
- SOPs for compensation settings

A red starburst graphic with multiple points, containing white text.

**More objective,  
reproducible,  
high quality in  
results between  
different labs**

- Preselected - previous experience:

FITC, PE, PerCP, PerCP-Cy5.5, PECy7

APC


- Extensive testing - based on

- ✓ MFI
- ✓ Stain index
- ✓ Spillover
- ✓ Stability
- ✓ Availability

APC-Cy7, AF700, APC-H7,

PacB, HV450, PacO, AmCyan, HV500

# EuroFlow SOP for sample preparation



Add antibodies

- Multiple tubes? First Abs



Wash buffer  
– final volume 50  $\mu$ L

- **NaN<sub>3</sub> – Sodium azide**
  - Inhibits potential redistribution of surface antigens
- **BSA – Bovine serum albumin**
  - Prevents non-specific binding

Add 50  $\mu$ L of sample

- Final volume 100  $\mu$ L
- **Incubate 30 min**

FACS Lysing solution

- 1:10 distilled water
- **Incubate 10 min**

Spin → **wash** → spin

- 540g, 5 min
- **Resuspend** 2mL (no NaN<sub>3</sub>)

Resuspend → **Acquire**

- Within an hour! 4°C, dark





# EuroFlow SOPs: combined intracellular and surface membrane staining

- **ALOT, MM MRD tube 2, PCST**

## Membrane staining



Fix

Wash → spin

Permeabilize

Stain

Add antibodies

- Multiple tubes? First Abs



Wash buffer  
– final volume 50 µL

- **NaN<sub>3</sub> – Sodium azide**
  - Inhibits potential redistribution of surface antigens
- **BSA – Bovine serum albumin**
  - Prevents non-specific binding

Add 50 µL of sample

- Final volume 100 µL
- **Incubate 30 min**


- Solution A, **MIX**
- 15 min, RT, dark



- 540g, 5 min
- Resuspend

- Solution B, **MIX**
- 15 min, RT, dark

- Intracellular
- 15 min, RT, dark



Spin -> **wash**-> spin

- 540g, 5 min
- **Resuspend** 2mL wash buffer

Resuspend -> **Acquire**


- Without NaN<sub>3</sub>
- Within an hour! 4°C, dark

# EuroFlow SOPs – BulkLysis protocol



Why? To increase the sensitivity of MRD detection or when higher nucleated cell concentration is needed

- ✓ MRD in Multiple Myeloma (MM)
- ✓ MRD in B-cell Precursor Acute Lymphoblastic Leukemia (BCP-ALL)
- ✓ Primary Immunodeficiencies Orientation Tube (PIDOT)



# EuroFlow SOPs – BulkLysis protocol

Sample in 50 mL tube


- Max 2mL in one tube
- MRD: min  $10 \times 10^6$  nucleated cells per tube



Fill the tube up to 50 mL with Bulk lyse

- **diluted** to 1X in dH2O at RT

15 min roller



800g, 10min → discard supernatant


- ~ 300  $\mu$ L remains

Resuspend → Wash

- Add 2 mL wash buffer → **resuspend** → fill up to 50 mL wash buffer

800g, 5min → discard supernatant

- ~ 300  $\mu$ L remains





# EuroFlow SOPs – BulkLysis protocol

Resuspend → Transfer to FACS tube

- Add 2 mL wash buffer → **resuspend**




Recover the rest of the cells

- Wash the Falcon tube with 2 mL wash buffer → transfer to the FACS tube



540g, 5 min → discard supernatant

- If more tubes were used → combine
- Final volume  $\geq 300 \mu\text{L}$



Adjust concentration!

- **$1 \times 10^5$  cells/ $\mu\text{L}$**



## FOLLOW APPROPRIATE **STAINING PROTOCOL**:

- PIDOT, BCP-ALL MRD, MM MRD Tube 1: membrane staining
- MM MRD tube 2: intracellular staining



## SOP for standardized instrument settings – FSC + SSC

BD FACSCanto II:

FSC: 55,000 (range 50,000 – 60,000)

SSC: 13,000 (range 11,000 – 15,000)


Threshold in FSC: 10,000






- Gating
- Comparison
- Follow up
- Database
- Merge

5x normal PB stained with LST, merged

## A. Record time parameter to detect technical issues



System fluctuations (laser, sheath pressure)?  
Cell aggregation/clumping?  
Bubbles?



## E. Anticoagulants affect the expression of certain markers





- EuroFlow method has been validated with samples collected with **EDTA** anticoagulant (purple top tube)
- The use of other anticoagulants, e.g. **heparin**, affects the expression of certain markers



Downregulation of CD138 expression in the membrane of plasma cells after heparin exposure

*Flores-Montero J, de Tute R, Paiva B, Perez JJ, Böttcher S, Wind H, Sanoja L, Puig N, Lecrevisse Q, Vidriales MB, van Dongen JJM and Orfao A. Immunophenotype of Normal vs. Myeloma Plasma Cells: Toward Antibody Panel Specifications for MRD Detection in Multiple Myeloma. Cytometry Part B 2016; 90B: 61–72.*

## D. Compensation issues



Use single-stained controls (cells or beads) with **same fluorochrome** as used in the assay  
Unstained control must have **same autofluorescence** as single-stained controls  
Tandem dyes need **lot specific compensation**

# Plasmatic Cells Screening Tube



- **Consensus 8-color panel:**

PacB	OC515	FITC	PE	PerCP-Cyanine5.5	PE-Cyanine7	APC	APC-C750
CD45	CD138	CD38	CD56	$\beta_2$ microglobulin	CD19	CylgKappa	CylgLambda
CD45*	CD138*	CD38	CD28	CD27	CD19	CD117	CD81

V450	V500c	FITC	PE	PerCP-Cyanine5.5	PE-Cyanine7	APC	APC-C750
CD138	CD45	CD38	CD56	$\beta_2$ microglobulin	CD19	CylgKappa	CylgLambda
CD138	CD45	CD38	CD28	CD27	CD19	CD117	CD81

- Identification, discrimination and enumeration of normal and aberrant plasma cells
- Initial evaluation and screening of PCD in BM samples
- $B_2$ -microglobulin allows for prognostic stratification



**ORIGINAL ARTICLE**

# Next Generation Flow for highly sensitive and standardized detection of minimal residual disease in multiple myeloma

J Flores-Montero<sup>1,19</sup>, L Sanoja-Flores<sup>1,19</sup>, B Paiva<sup>2,19</sup>, N Puig<sup>3</sup>, O García-Sánchez<sup>3</sup>, S Böttcher<sup>4</sup>, VHJ van der Velden<sup>5</sup>, J-J Pérez-Morán<sup>3</sup>, M-B Vidriales<sup>3</sup>, R García-Sanz<sup>3</sup>, C Jimenez<sup>3</sup>, M González<sup>3</sup>, J Martínez-López<sup>6</sup>, A Corral-Mateos<sup>1</sup>, G-E Grigore<sup>7</sup>, R Fluxá<sup>7</sup>, R Pontes<sup>8</sup>, J Caetano<sup>9</sup>, L Sedek<sup>10</sup>, M-C del Cañizo<sup>3</sup>, J Bladé<sup>11</sup>, J-J Lahuerta<sup>6</sup>, C Aguilar<sup>12</sup>, A Báez<sup>13</sup>, A García-Mateo<sup>14</sup>, J Labrador<sup>15</sup>, P Leoz<sup>1</sup>, C Aguilera-Sanz<sup>16</sup>, J San-Miguel<sup>2,20</sup>, M-V Mateos<sup>3,20</sup>, B Durie<sup>17,21</sup>, JJM van Dongen<sup>5,18,21</sup> and A Orfao<sup>1,21</sup>

# Reagent kit for MM MRD/CTPC



- **Consensus 8-color 2-tube panel:**


MM MRD	BV421	BV510	FITC	PE	PerCP-Cyanine 5.5	PE-Cyanine7	APC	APC-C750
Tube 1	CD138*	CD27*	CD38-multiepitope	CD56	CD45	CD19	CD117	CD81
Tube 2	CD138*	CD27*	CD38-multiepitope	CD56	CD45	CD19	CyIgKappa	CyIgLambda

- High sensitivity of close to  $10^{-6}$  as recommended by the IMWG
- CD38 multiepitope to overcome antigen loss post-drug treatment
- Replicate consistency by the use of two tubes
- Pre-mix lyophilized format to avoid pipetting errors and to reduce time and costs
- Ancillary reagents included: BulkLysis™, FIX&PERM® , specific compensation vials
- Minimally-invasive test to evaluate tumor load in peripheral blood

*\*CD138-BV421 and CD27-BV510 to be added as drop-ins*





## H. MM MRD: different processing Tube 1 vs tube 2



Median values from both tubes should look a bit different


# Treatment considerations: CD38 Multiepitope

Targets CD38 molecule on plasma cells, even in patients under anti-CD38 therapy.




# Número de células colectadas es CRÍTICO para la sensibilidad del ensayo


CD45 vs SSC-A




100 000 events




500 000 events



1 million events








5 million events



10 million events

CD45 vs CD38





# Sensitivity


IMWG recommendation: Sensitivity  $10^{-5}$




$$\% LOD = 100 \times \frac{20}{\text{Nucleated cells}}$$

$$\% LOQ = 100 \times \frac{50}{\text{Nucleated cells}}$$

$$\% cPC = 100 \times \frac{cPC}{\text{Nucleated cells}}$$



% LOD	0.0002	0.0004	0.0008
% LOQ	0.0005	0.0010	0.0021
% cPC	<b>0.0031</b>	<b>0.0016</b>	<b>0.0008</b>






% LOD	0.0021	0.0041	0.0207
% LOQ	0.0052	0.0104	0.0518
% cPC	<b>0.0003</b>	<b>0.0002</b>	<b>0.0000</b>

# 2016: MRD monitoring using 2<sup>nd</sup> generation flow affords higher sensitivity


## 1<sup>st</sup> generation flow

2X10<sup>5</sup> cells analyzed



## 2<sup>nd</sup> generation flow

2X10<sup>6</sup> cells analyzed




1<sup>st</sup> generation flow

## 2<sup>nd</sup> generation flow



	-	+
-		15/50 (30%)
+	0/50 (0%)	35/50 (70%)

The gating strategy for the identification of a total plasma cell population and the discrimination of normal vs clonal plasma cells. An example of the sequential steps followed to gate the total plasma cell population are shown in Panel A. Panel B shows the discrimination of the clonal plasma cell population in this sample based on their abnormal phenotypic profile (CD38<sup>lo</sup>, CD45<sup>-</sup>, CD19<sup>-</sup>, CD56<sup>het+</sup>, CD81<sup>-</sup>, CD117<sup>+</sup> and  $\text{cyIg}^{\lambda^+} / \text{cyIg}^{\kappa^-}$ ). The total plasma cell population corresponds to 0.55% of BM cellularity while the clonal plasma cells represent 0.27%. Please note that the phenotypic profile of clonal plasma cells might be different in each patient sample.





Next gen flow-minimal residual disease in myeloma  
J Flores-Montero *et al*



Parameter	% of contribution
CD56	28%
CD19	19%
CD45	17%
CD81	16%
CD138	6%
CD27	4%

# AUTOMATIC REPORT INFORMATION

## Cellularity

Frequency of the populations in the bone marrow of study (% of events within the reference population). In bold, populations altered compared to the reference range.

Population	Frequency	Reference
<b>Plasma cells</b>	<b>0.029</b>	(0.048 - 0.97)
<b>B cells</b>	<b>0.38</b>	(1 - 5.2)
<b>Mature B cells</b>	<b>0.38</b>	(0.95 - 3.6)
<b>T and NK and basophils</b>	<b>25.7</b>	(8.8 - 17.4)
Eosinophils	1.1	(0.68 - 2.3)

## Other features:

**Absent populations:** populations not present in the sample of study at the moment the report is opened. They can be in the unchecked events or not present at all. These populations should be found in a normal sample and are included in the database.

Although the analysis displays all absent populations, the report informs only if they are B-cell precursors, Mast cells, B cells, Mature B cells, T and NK and Basophils, Eosinophils, Neutrophils, Monocytes, Myeloid precursors, Nucleated red cells.

Absent Populations		7	0
<input type="checkbox"/> CD19-CD56+ PC		3	0
<input type="checkbox"/> CD19-CD56- PC		3	0
<input type="checkbox"/> B-cell precursors		1	0




## LYMPHOID NEOPLASIA

# Detailed characterization of multiple myeloma circulating tumor cells shows unique phenotypic, cytogenetic, functional, and circadian distribution profile

Bruno Paiva,<sup>1</sup> Teresa Paino,<sup>2</sup> Jose-Maria Sayagues,<sup>2,3</sup> Mercedes Garayoa,<sup>2</sup> Laura San-Segundo,<sup>2</sup> Montserrat Martín,<sup>2</sup> Ines Mota,<sup>3</sup> María-Luz Sanchez,<sup>2,3</sup> Paloma Bárcena,<sup>2,3</sup> Irene Aires-Mejia,<sup>2</sup> Luis Corchete,<sup>2</sup> Cristina Jimenez,<sup>2</sup> Ramon Garcia-Sanz,<sup>2</sup> Norma C. Gutierrez,<sup>2</sup> Enrique M. Ocio,<sup>2</sup> Maria-Victoria Mateos,<sup>2</sup> Maria-Belen Vidriales,<sup>2</sup> Alberto Orfao,<sup>2,3</sup> and Jesús F. San Miguel<sup>1</sup>


<sup>1</sup>Clinica Universidad de Navarra, Centro de Investigacion Medica Aplicada, Pamplona, Spain; <sup>2</sup>Hospital Universitario de Salamanca, Instituto de Investigaion Biomedica de Salamanca, Instituto de Biologia Molecular y Celular del Cancer (Univesidad de Salamanca - Consejo Superior de Investigaciones Cientificas), Salamanca, Spain; and <sup>3</sup>Servicio General de Citometría and Departamento de Medicina, Universidad de Salamanca, Salamanca, Spain



CSC, cancer stem cell; IGH, immunoglobulin heavy chain; MGUS, monoclonal gammopathy of undetermined significance; MM, multiple myeloma.


1. Kyle RA, et al. N Engl J Med. 2002;346:564-569.
2. Kyle RA, et al. N Engl J Med. 2007;356:2582-2590.

# Circulating tumor cells (CTCs) progressively increase from MGUS to MM




# Clinical significance of CTCs quantification in MGUS and smoldering MM

MGUS <sup>1</sup>




Smoldering MM <sup>2</sup>




1. Kumar. J Clin Oncol. 2005 20;23(24):5668-74
2. Bianchi G, et al. Leukemia. 2013;27(3):680-5

# CTCs and CD34+ HSCs share the same circadian rhythm regulated by the SDF1/CXCR4 axis

MM-CTCs (median cells/ $\mu$ L)  
CD34+ HSC (median cells/ $\mu$ L)




CXCR4 (Amount of antigen MFI expression / MM-CTC)  
SDF-1 $\alpha$  levels (pg/mL)



MM patients at relapse (n=6)  
Quantification started at 16:00pm every 4h up to 12:00am next day (when patients' initiated treatment)  
Time points 16h and 21h have been duplicated to facilitate viewing of the time curve


# Multiple myeloma: a phenotypic perspective






# The potential to egress into PB is restricted to a minor subclone in the BM...

BM MM-PC vs. CTCs: principle component analysis (APS) of 22 antigens



*...with an unique profile of integrin and adhesion molecules*

# Diagnosis and Monitoring by Flow Cytometry





## Diagnosis

- Disease staging
- Cell level light chain restriction

## Prognosis

- Identification of patients with sustained disease control due to high normal PC recovery and an immune profile with increased B-cell maturation
- Implementation of individualized therapy monitoring strategies


## MRD Assessment

- Detection of abnormal plasma cells at very low levels by highly sensitive technologies
- Compare the efficacy of different treatment strategies

# NEXT GENERATION FLOW (NGF)

## Clinical Question

<b>Panel Design</b>	Meaningful markers Best Conjugates
<b>Instrument Setup</b>	Calibration Compensation
<b>Sample Processing</b>	BulkLysis™ for sample concentration Standard operating procedures (SOP)
<b>Data Acquisition</b>	Required sensitivity - number of events acquired Acquisition speed
<b>Data Analysis</b>	Precise identification of all normal cells Characterization of abnormal phenotypes and distributions



# ADVANTAGES NGF

## Flow cytometry evaluation of PCDs

Characterization of the **immune cell distribution** and function.


Orientation to specific genetic tests and obtaining of a “probable diagnosis”.

- ✓ Efficient.
- ✓ Widely accessible and minimally invasive.
- ✓ Fast (approximately 4h TAT).
- ✓ Minimizes the need of complete clinical information.

## NGF approach for MM MRD evaluation

- ✓ Highly sensitive (20 cells in  $10^6$ ) even for those populations with lower frequency.
- ✓ Objective (automated analysis and reporting).
- ✓ Extensively validated by EuroFlow™.
- ✓ External QC assessment available for laboratory performance and sample processing.

# EUROFLOW PANELS FOR HEMATO-ONCOLOGY



# INFINCYT™: ANALYSIS SOFTWARE











# MM-MRD Manual Analysis

## MANUAL ANALYSIS

- X Subject to Interpretation
- X Multidimensional Complexity
- X Time Consuming




# MM MRD automated analysis in Infinicyt


*AUTOMATED  
GATING & IDENTIFICATION*

*EXPERT USER*



*AUTOMATIC REPORT*



Raw MM MRD flow data



Identified cell populations



Infinicyt identifies normal populations in the sample

The user reviews and confirms what is not normal  
(e.g. aberrant plasma cells)

## CELLULARITY (Data referred to 100% of the viable cells)

Population	Frequency	Reference	Population	Frequency	Reference
Plasma cells	0.094	(0.048 - 0.97)	Neutrophils	77.3	(60.1 - 78.6)
B cells	1.4	(1 - 5.2)	Monocytes	4.6	(3.4 - 8.1)
Mature B cells	1	(0.95 - 3.6)	Myeloid precursors	1.1	(0.92 - 3.4)
T and NK and basophils	7.7	(8.8 - 17.4)	Mast cells	0.0011	(0.0048 - 0.024)
Eosinophils	1.4	(0.68 - 2.3)	Nucleated red cells	6.8	(2.2 - 15.8)
			Abnormal plasma cells	0.017	-

## Resultados estadísticos con referencia de poblaciones normales

Limit of Detection [LOD]:	0.00024	Lower Limit of Quantitation [LLQ]:	0.00061
---------------------------	---------	------------------------------------	---------

sensibilidad

nuevos conceptos!!

LOD detectable 20

LOQ cuantificable 50

descripción  
inmunofenotípica  
de las CPA

### ABNORMAL PLASMA CELLS IMMUNOPHENOTYPE

FSC<sup>normal</sup> het SSC<sup>normal</sup> CD138<sup>+</sup>  
CD38<sup>lo</sup> het CD56<sup>+</sup> CD45<sup>++</sup> CD19<sup>-</sup> CD27<sup>lo</sup> het CD117<sup>+</sup> CD81<sup>+</sup>  
CylgLambda<sup>+</sup>

- marcadores normales de las CPA
- marcadores alterados
- clonalidad

### COMMENT

Plasma cells are present at 0.11% of total nucleated cells; 15.3% of them express an aberrant phenotype (CD38<sup>lo</sup> het CD56<sup>+</sup> CD45<sup>++</sup> CD19<sup>-</sup> CD27<sup>lo</sup> het CD117<sup>+</sup> CD81<sup>+</sup>).  
Sample probably hemodiluted (there is a decreased number of mast cells).  
Normal plasma cells have altered Cylg K/L ratio. Confirm/discard the normality of these cells.

comentarios  
importantes

### CONCLUSION

BM compatible with POSITIVE MRD. Note: An hemodiluted sample will underestimate or even fail to detect the percentage of BM infiltration.

## AUTOMATED GATING & IDENTIFICATION REPORT

# MM MRD automated analysis in Infinicyt

### CELLULARITY (Data refers to 100% of viable cells)

Population	Frequency	Reference	Population	Frequency	Reference
<b>Plasma cells (PC)</b>	<b>0.019</b>	(0.048 - 0.97)	<b>Neutrophils</b>	<b>58.6</b>	(60.1 - 78.6)
B cells	2.7	(1 - 5.2)	<b>Monocytes</b>	<b>8.4</b>	(3.4 - 8.1)
Mature B cells	1.2	(0.95 - 3.6)	Myeloid precursors ^	1.4	(0.92 - 3.4)
B-cell precursors	1.5	(0.085 - 1.9)	<b>Mast cells ^</b>	<b>0.004</b>	(0.0048 - 0.024)
<b>T and NK and basophils</b>	<b>20.6</b>	(8.8 - 17.4)	Nucleated red cells	7	(2.2 - 15.8)
Eosinophils	2	(0.68 - 2.3)	<b>Abnormal PC</b>	<b>0.0073</b>	-

^Frequency calculation based on a single file.

Limit of Detection (LOD):	0.00022	Lower Limit of Quantification (LLOQ):	0.00056
---------------------------	---------	---------------------------------------	---------

### IMMUNOPHENOTYPE OF ABNORMAL PLASMA CELLS

CD81<sup>+</sup> CD138<sup>+</sup>  
 CD117<sup>+</sup> CD27<sup>+</sup> CD38<sup>lo</sup> FSC<sup>hi</sup> <sup>het</sup>CD19<sup>+</sup> CD56<sup>+</sup> CD45<sup>SSC</sup><sup>hi</sup>  
 CylgKappa<sup>+</sup>

### COMMENT

0.026% of Plasma cells were detected: 27.9% of these express an aberrant phenotype ( CD117<sup>+</sup> CD27<sup>+</sup> CD38<sup>lo</sup> FSC<sup>hi</sup> <sup>het</sup>CD19<sup>+</sup> CD56<sup>+</sup> CD45<sup>SSC</sup><sup>hi</sup>).  
 Sample possibly contaminated with peripheral blood (decreased percentage of mast cells was observed).

### CONCLUSION

Bone marrow compatible with positive MRD. Note: Possible hemodilution of the sample must be considered to assess the percentage of infiltration.

### MM-MRD Assessment

- Abnormal Plasma Cells

### Quality of the Sample: Hemodilution

- Mast Cells


### Sensitivity

- Number of events

### Bone Marrow Immune Regeneration

- B-Cell Precursors
- Nucleated Red Cells

# NGF reaches $10^{-6}$ sensitivity in the vast majority of MM patients (GEM2012MENOS65)



LOQ (level of quantitation) : 50 cells / total nucleated viable cells

LOD (level of detection): 20 cells / total nucleated viable cells


Flores-Montero J, et al. Leukemia. 2017. doi: 10.1038/leu.2017.29

**LOQ: límite de cuantificación: 50 células/ total de células nucleadas viables**


**LOD: límite de detección: 20 células/ total de células nucleadas viables**



# PATIENTS IN VGPR, CR and sCR




# PATIENTS IN CR and sCR



# 2016: MRD monitoring using 2<sup>nd</sup> generation flow improves discrimination between normal vs. clonal PCs


## 4-color flow-MRD




82% accuracy  
(41/50 patients)

<input checked="" type="checkbox"/>	CD38	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	CD138	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	CD19	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	CD56	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	CD27	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	CD81	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	CD117	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	CD45	<input checked="" type="checkbox"/>

## 8-color flow-MRD





96% accuracy  
(48/50 patients)



- Muestra óptima: médula ósea con EDTA, sin hemodilución
- Cantidad de células disponibles: adquirir 10 millones!!!

# Next generation Flow-MRD monitoring in MM

*- applicable to all patients -*





## **Conclusions**

**-MRD positive or negative:** Compatible with MM MRD positive conclusion will be reported for values above the LOD.

**-Warning about the quality of the sample (hemodilution):** The criteria used for hemodilution is the frequency of mast cells lower than the range of the database.

Note: Mesenchymal cells were also considered, but finally were discarded because this panel could not accurately identify them. Nucleated red cells were also considered but finally discarded because there is a great heterogeneity within normal samples. Mast cells population is more stable and homogeneous, it can be affected by treatment but it is not very common.

**-Warning about the number of events acquired:** Recommended standard about number of events to be acquired is 10 Million to reach a sensitivity of  $10^{-5}$ . It takes into account an acquisition of 10 Millions with typically 20% of debris. The warning appears if less than 8 Million of nucleated cells if found at the end of analysis.



**Gracias por su atención!!!**